

Simultaneous Analysis of Volatile Organic Compounds in Water Using Purge and Trap Gas-Chromatography–Mass-Spectrometry

Tatsuya Shimofuji, Kazuyuki Ishii, Noriyuki Baba, Hiroshi Hayashida,
GL Sciences Inc. Tokyo, Japan

Introduction

The analysis conditions for purge and trap (P&T) gas-chromatograph–mass-spectrometer (GC-MS) are specified in US EPA Method 524.2 for volatile organic compounds (VOC) analysis in drinking water. A trap tube packed with carbon is used in this method. As a result, reducing the amount of water introduced into the GC-MS during desorption is effective for performing stable analysis. In this application note, we report on the simultaneous analysis of 83 VOC components using the PT7000 system, which is a purge and trap concentration introduction system equipped with a high-performance MCS.

Method

The purge and trap system included PT7000 and vial autosampler AS7100, which have been used for highly sensitive analysis of mold odor substances. GC-MS used the GCMS-QP2020 NX (Shimadzu Corporation). The MCS of the PT7000 was subjected to a highly inert treatment to reduce the influence of VOC adsorption, and achieve efficient and stable water removal. For VOC analysis, the GC capillary column used was a moderately polar column TC-VMS (0.25 mm I.D. × 30 m df = 1.4 μm). The desorb time was reduced from 4 to 2 min, as recommended by EPA Method 524.2. Tables 1 and 2 show the detailed analytical conditions, including other parameters. In addition to Restek's Calibration Mix # 4, # 5, # 7A, # 7B, and # 8A, the VOC standard solution was made by dissolving three components of Tokyo Kasei Industry's bromomethane, dichlorodifluoromethane, and trichlorofluoromethane in methanol, mixing them, and preparing 100 ppm as a mixed standard stock solution of 83 components. The calibration curve standard solutions were prepared in 40 mL vials by adding the mixed standard stock solution to the water sample in increments of 0.5 to 10 ppb, respectively. Fluorobenzene was used as an internal standard, and *p*-bromofluorobenzene and 1,2-dichlorobenzene-*d*₄ were used as surrogates. They were mixed so that their concentrations after addition in water were 5 ppb, respectively, and they were placed in an internal bottle, and 2 μL per 5 mL of the sample was automatically added using the vial autosampler function.



Table 3 Calibration, Accuracy and Precision

| No. | Compounds | Calibration (0.1 - 10 ppb) | | MDL (n=7, 0.5 ppb) | | | |
|-------|--|----------------------------|-----------------------------|--------------------|-----------|------------------|-------------------|
| | | R.T. (min) | Linearity (r ²) | Ave. Conc. (ppb) | MDL (ppb) | Accuracy (± 20%) | Precision (≤ 20%) |
| 1 | Dichlorodifluoromethane | 0.919 | 0.9993 | 0.47 | 0.123 | 94 | 8.7 |
| 2 | Chloromethane | 1.031 | 0.9985 | 0.46 | 0.103 | 91 | 7.5 |
| 3 | Vinyl chloride | 1.063 | 0.9997 | 0.48 | 0.138 | 95 | 9.7 |
| 4 | Bromomethane | 1.234 | 0.9999 | 0.41 | 0.211 | 82 | 17.2 |
| 5 | Chloroethane | 1.302 | 0.9997 | 0.46 | 0.130 | 92 | 9.4 |
| 6 | Trichlorofluoromethane | 1.376 | 1.0000 | 0.48 | 0.125 | 96 | 8.7 |
| 7 | Ethyl ether | 1.566 | 1.0000 | 0.47 | 0.119 | 94 | 8.5 |
| 8 | 1,1-Dichloroethylene | 1.671 | 1.0000 | 0.47 | 0.134 | 94 | 9.5 |
| 9 | Carbon disulfide | 1.683 | 0.9999 | 0.48 | 0.066 | 96 | 4.6 |
| 10 | Iodomethane | 1.755 | 0.9991 | 0.40 | 0.242 | 79 | 20.4 |
| 11 | Allyl chloride | 1.982 | 0.9999 | 0.47 | 0.071 | 95 | 5.0 |
| 12 | Dichloromethane | 2.055 | 0.9997 | 0.48 | 0.089 | 96 | 6.2 |
| 13 | Acetone | 2.115 | 0.9981 | 0.50 | 0.192 | 101 | 12.7 |
| 14 | <i>trans</i> -1,2-Dichloroethylene | 2.167 | 0.9999 | 0.47 | 0.075 | 95 | 5.3 |
| 15 | MTBE | 2.283 | 0.9999 | 0.47 | 0.063 | 94 | 4.4 |
| 16 | 1,1-Dichloroethane | 2.671 | 0.9999 | 0.47 | 0.097 | 94 | 6.9 |
| 17 | Acrylonitrile | 2.738 | 0.9994 | 0.46 | 0.153 | 93 | 11.0 |
| 18 | <i>cis</i> -1,2-Dichloroethylene | 3.165 | 0.9999 | 0.49 | 0.049 | 98 | 3.3 |
| 19 | 2,2-Dichloropropane | 3.266 | 0.9998 | 0.50 | 0.030 | 101 | 2.0 |
| 20 | Bromochloromethane | 3.351 | 1.0000 | 0.49 | 0.050 | 99 | 3.3 |
| 21 | Chloroform | 3.448 | 0.9998 | 0.46 | 0.023 | 93 | 1.6 |
| 22 | Tetrachloromethane | 3.551 | 0.9998 | 0.48 | 0.067 | 96 | 4.7 |
| 23 | Tetrahydrofuran | 3.635 | 0.9994 | 0.46 | 0.025 | 92 | 1.8 |
| 24 | 2-Butanone | 3.634 | 0.9993 | 0.54 | 0.047 | 108 | 2.9 |
| 25 | 1,1,1-Trichloroethane | 3.631 | 0.9999 | 0.48 | 0.045 | 96 | 3.1 |
| 26 | Methyl acrylate | 3.635 | 0.9989 | 0.48 | 0.034 | 95 | 2.4 |
| 27 | 1,1-Dichloropropene | 3.767 | 0.9997 | 0.48 | 0.026 | 95 | 1.8 |
| 28 | 1-Chlorobutane | 3.839 | 0.9994 | 0.48 | 0.025 | 96 | 1.7 |
| 29 | Benzene | 4.034 | 0.9998 | 0.48 | 0.024 | 95 | 1.6 |
| 30 | Propionitrile | 4.115 | 0.9999 | 0.50 | 0.052 | 99 | 3.5 |
| 31 | Methacrylonitrile | 4.132 | 0.9992 | 0.46 | 0.056 | 92 | 4.1 |
| 32 | 1,2-Dichloroethane | 4.273 | 0.9999 | 0.48 | 0.033 | 97 | 2.3 |
| I.S. | Fluorobenzene (I.S.) | 4.527 | - | - | - | - | - |
| 33 | Trichloroethene | 4.722 | 0.9999 | 0.48 | 0.039 | 96 | 2.7 |
| 34 | Dibromomethane | 5.230 | 1.0000 | 0.49 | 0.049 | 98 | 3.3 |
| 35 | 1,2-Dichloropropane | 5.377 | 0.9998 | 0.50 | 0.055 | 100 | 3.7 |
| 36 | Bromodichloromethane | 5.490 | 0.9998 | 0.49 | 0.048 | 97 | 3.3 |
| 37 | Methyl methacrylate | 5.820 | 0.9981 | 0.46 | 0.076 | 92 | 5.5 |
| 38 | <i>cis</i> -1,3-Dichloropropene | 6.383 | 0.9995 | 0.48 | 0.040 | 96 | 2.8 |
| 39 | Toluene | 6.720 | 0.9993 | 0.47 | 0.042 | 94 | 3.0 |
| 40 | Chloroacetonitrile | 7.010 | 0.9851 | 0.53 | 0.147 | 105 | 9.3 |
| 41 | 2-Nitropropane | 7.069 | 0.9995 | 0.49 | 0.040 | 97 | 2.8 |
| 42 | 1,1-Dichloropropanone | 7.078 | 0.9997 | 0.46 | 0.060 | 92 | 4.3 |
| 43 | Tetrachloroethylene | 7.220 | 0.9998 | 0.49 | 0.028 | 98 | 1.9 |
| 44 | <i>trans</i> -1,3-Dichloropropene | 7.325 | 0.9992 | 0.47 | 0.036 | 94 | 2.5 |
| 45 | 4-Methyl-2-pentanone | 7.330 | 0.9988 | 0.44 | 0.033 | 89 | 2.5 |
| 46 | 1,1,2-Trichloroethane | 7.499 | 1.0000 | 0.49 | 0.056 | 97 | 3.8 |
| 47 | Ethyl methacrylate | 7.616 | 0.9975 | 0.44 | 0.070 | 87 | 3.8 |
| 48 | Dibromochloromethane | 7.671 | 0.9996 | 0.48 | 0.051 | 95 | 3.6 |
| 49 | 1,3-Dichloropropane | 7.790 | 0.9998 | 0.47 | 0.037 | 94 | 2.6 |
| 50 | 1,2-Dibromoethane | 7.887 | 0.9998 | 0.47 | 0.037 | 95 | 2.6 |
| 51 | 2-Hexanone | 8.259 | 0.9968 | 0.44 | 0.136 | 89 | 10.2 |
| 52 | Chlorobenzene | 8.456 | 0.9998 | 0.47 | 0.048 | 95 | 3.4 |
| 53 | Ethylbenzene | 8.524 | 0.9990 | 0.45 | 0.052 | 91 | 3.8 |
| 54 | 1,1,1,2-Tetrachloroethane | 8.538 | 0.9999 | 0.49 | 0.062 | 97 | 4.3 |
| 55 | <i>m,p</i> -Xylene | 8.672 | 0.9992 | 0.45 | 0.054 | 90 | 4.0 |
| 56 | <i>o</i> -Xylene | 9.052 | 0.9991 | 0.47 | 0.046 | 94 | 3.3 |
| 57 | Bromoform | 9.084 | 0.9993 | 0.49 | 0.036 | 99 | 2.4 |
| 58 | Styrene | 9.101 | 0.9979 | 0.42 | 0.058 | 84 | 4.6 |
| 59 | Isopropylbenzene | 9.339 | 0.9986 | 0.46 | 0.086 | 92 | 6.2 |
| SURR. | <i>p</i> -Bromofluorobenzene (SURR.) | 9.551 | - | - | - | - | - |
| 60 | Bromobenzene | 9.619 | 0.9972 | 0.50 | 0.058 | 103 | 3.8 |
| 61 | <i>n</i> -Propylbenzene | 9.687 | 0.9986 | 0.50 | 0.084 | 101 | 5.6 |
| 62 | 1,1,2,2-Tetrachloroethane | 9.752 | 0.9953 | 0.51 | 0.043 | 103 | 2.8 |
| 63 | 2-Chlorotoluene | 9.787 | 0.9979 | 0.51 | 0.045 | 102 | 2.9 |
| 64 | 1,2,3-Trichloropropane | 9.840 | 0.9957 | 0.52 | 0.095 | 105 | 6.0 |
| 65 | 1,3,5-Trimethylbenzene | 9.859 | 0.9996 | 0.50 | 0.107 | 100 | 7.2 |
| 66 | <i>trans</i> -1,4-Dichloro-2-butene | 9.893 | 0.9990 | 0.48 | 0.053 | 96 | 3.6 |
| 67 | 4-Chlorotoluene | 9.926 | 0.9965 | 0.49 | 0.089 | 98 | 6.0 |
| 68 | <i>tert</i> -Butylbenzene | 10.103 | 0.9978 | 0.51 | 0.082 | 101 | 5.4 |
| 69 | Pentachloroethane | 10.100 | 0.9928 | 0.58 | 0.188 | 116 | 10.8 |
| 70 | 1,2,4-Trimethylbenzene | 10.161 | 0.9992 | 0.49 | 0.103 | 99 | 7.0 |
| 71 | <i>sec</i> -Butylbenzene | 10.244 | 0.9990 | 0.50 | 0.083 | 99 | 5.6 |
| 72 | <i>p</i> -Isopropyltoluene | 10.365 | 0.9993 | 0.49 | 0.108 | 99 | 7.3 |
| 73 | <i>m</i> -Dichlorobenzene | 10.383 | 0.9971 | 0.51 | 0.063 | 103 | 4.1 |
| 74 | <i>p</i> -Dichlorobenzene | 10.457 | 0.9981 | 0.51 | 0.061 | 101 | 4.0 |
| 75 | <i>n</i> -Butylbenzene | 10.686 | 0.9991 | 0.49 | 0.070 | 99 | 4.7 |
| 76 | Hexachloroethane | 10.753 | 0.9978 | 0.52 | 0.138 | 103 | 8.9 |
| SURR. | 1,2-Dichlorobenzene- <i>d</i> ₄ (SURR.) | 10.767 | - | - | - | - | - |
| 77 | 1,2-Dichlorobenzene | 10.774 | 0.9984 | 0.51 | 0.076 | 101 | 5.0 |
| 78 | 1,2-Dibromo-3-Chloropropane | 11.381 | 0.9907 | 0.51 | 0.155 | 102 | 10.1 |
| 79 | Nitrobenzene | 11.797 | 0.9949 | 0.55 | 0.332 | 111 | 20.0 |
| 80 | Hexachlorobutadiene | 11.863 | 0.9957 | 0.53 | 0.047 | 107 | 3.0 |
| 81 | 1,2,4-Trichlorobenzene | 11.880 | 0.9997 | 0.53 | 0.119 | 107 | 7.5 |
| 82 | Naphthalene | 12.120 | 0.9988 | 0.51 | 0.077 | 103 | 5.0 |
| 83 | 1,2,3-Trichlorobenzene | 12.250 | 0.9989 | 0.49 | 0.034 | 98 | 2.4 |

Table 1 P&T Method Conditions

| Standby | Bake |
|----------------------|-----------------|
| Valve oven temp. | 150 °C |
| Transfer line temp. | 150 °C |
| Mount temp. | 60 °C |
| Purge ready temp. | 35 °C |
| MCS temp. | 40 °C |
| Standby flow rate | 40 mL/min |
| B.O.T temp. | 150 °C |
| Purge MCS temp. | 40 °C |
| Dry purge MCS temp. | 40 °C |
| Purge | |
| Prepurge time | 0 min |
| Prepurge flow rate | 0 mL/min |
| Sample temp. | 60 °C |
| Preheat time | 0.01 min |
| Purge time | 11 min |
| Purge temp. | 0 °C |
| Purge flow rate | 40 mL/min |
| Dry purge time | 1 min |
| Dry purge temp. | 20 °C |
| Dry purge flow rate | 40 mL/min |
| Desorb | |
| GC start | Start of desorb |
| Desorb Preheat temp. | 250 °C |
| Drain | ON |
| Desorb time | 2 min |
| Desorb temp. | 250 °C |
| Desorb flow rate | 40 mL/min |
| Bake time | 8 min |
| Bake temp. | 260 °C |
| MCS bake temp. | 260 °C |
| Bake flow rate | 80 mL/min |
| Bake sparger time | 10 min |
| After bake flow rate | 2 mL/min |
| After bake flow rate | 200 mL/min |
| Autosampler | |
| Sampling time | 0.2 min |
| I.S. addition | ON |
| Fill I.S. time | 0.2 min |
| Transfer time | 0.25 min |
| Loop wash time | 0.25 min |
| Loop purge time | 0.25 min |
| Rinse frequency | 3 |
| Bake fill time | 0.2 min |
| Bake transfer time | 0.25 min |
| Bake drain time | 0.25 min |
| Bake drain flow rate | 100 mL/min |
| Note | |
| Sparger | 5 mL |
| Trap | Vocarb 3000 |
| Sample loop | 5 mL |
| Chiller Tray | ON |
| Purge Gas | Nitrogen |

Table 2 GC - MS Method Conditions

| | |
|-------------|--|
| GC-MS | GCMS-QP2020 NX (Shimadzu Corporation) |
| Column | TC-VMS 0.25 mm I.D. × 30 m, df = 1.40 μm (GL Sciences Inc.) |
| Split ratio | 1:40 |
| Oven Temp. | 40 °C (hold 2.5 min) – 7 °C/min – 70 °C – 25 °C/min – 240 °C (1 min) |
| Carrier gas | He, 180 kPa |
| Scan range | 35 – 300 amu |




Fig.1 Total Ion Current Chromatogram (TICC)

Results and Conclusion

Fig. 1 depicts the total ion current chromatogram (TICC) obtained under the aforementioned analytical conditions. In 83 VOCs, excellent separation was demonstrated. The RSDs of the internal standard and surrogate were both less than 5%. Table 3 displays the retention time, calibration curve, and percent RSD at 0.5 ppb for each compound. As specified by the method, the percent RSD for all components were less than 20% (except Iodomethane). In fact, for most components, except 8 components, they were less than 10%. The linearity, MDL, and accuracy results were also satisfactory. The PT7000 system and GCMS-QP2020 NX combined analysis produced good results for VOC analysis in water and was applicable to EPA Method 524.

References

METHOD 524.2. MEASUREMENT OF PURGEABLE ORGANIC COMPOUNDS IN WATER BY CAPILLARY COLUMN GAS CHROMATOGRAPHY/MASS SPECTROMETRY
Revision 4.1, Tridited by J.W. Munch (1995)