

Notification of Shape Changes to MonoSpin S-type Spin Columns

Dear Valued Partner,

Thank you for your continuous support.

We are pleased to announce an update to our MonoSpin S-type spin columns, aimed at enhancing product quality.

The transition to the new shape will occur as current inventories are depleted. During this transition, please note that if you purchase multiple boxes, you may receive a mix of old and new versions in separate boxes. However, both versions will not be mixed within the same box.

Affected Products:

- MonoSpin S-type series 50/pk and 100/pk (Excluding Trypsin HP, ProA and ProG)
- MonoSpin Trialkit1,2,3 and Custom kit

For more details, please refer to the attached table at the end of this letter.

Details of Changes:





Change 1-1 Cap Closure Improvement: As indicated by the left arrow, the connection part has been modified to make cap closure easier.

Change 1-2 Diameter Adjustment: As indicated by the right arrow, the upper outer diameter has been reduced to prevent interference with adjacent spin columns during centrifugation.

Appearance Changes:

Please see the attached images for visual reference (Change 2).

The collection and waste tubes accompanying these products remain unchanged.

	Change 1	Change 2
【After Change】		
【Before Change】		

Impact on Quality and Reliability:

These modifications will not affect the quality or reliability of our products.

Implementation Timeline:

The switch to the new shape will take effect once current inventories are exhausted.

Thank you for your continued support of our products.

Yours sincerely,



Makoto Suzuki
 Manager, Quality Assurance Section
 GL Sciences Inc.

【Related Products】

Cat. No.	Descriptions
5010-21670	MonoSpin C18FF 50/pk
5010-21671	MonoSpin C18FF 100/pk
5010-21696	MonoSpin AG 50/pk
5010-21697	MonoSpin AG 100/pk
5010-21698	MonoSpin Phospholipid 50/pk
5010-21699	MonoSpin Phospholipid 100/pk
5010-21700	MonoSpin C18 50/pk
5010-21701	MonoSpin C18 100/pk
5010-21705	MonoSpin TiO 50/pk
5010-21706	MonoSpin TiO 100/pk
5010-21710	MonoSpin NH2 50/pk
5010-21711	MonoSpin NH2 100/pk
5010-21715	MonoSpin PBA 50/pk
5010-21716	MonoSpin PBA 100/pk
5010-21720	MonoSpin SAX 50/pk
5010-21721	MonoSpin SAX 100/pk
5010-21725	MonoSpin SCX 50/pk
5010-21726	MonoSpin SCX 100/pk
5010-21727	MonoSpin Amide 50/pk
5010-21728	MonoSpin Amide 100/pk
5010-21729	MonoSpin CBA 50/pk
5010-21730	MonoSpin CBA 100/pk
5010-21731	MonoSpin C18-CX 50/pk
5010-21732	MonoSpin C18-CX 100/pk
5010-21733	MonoSpin Ph 50/pk
5010-21734	MonoSpin Ph 100/pk
5010-21735	MonoSpin C18-AX 50/pk
5010-21736	MonoSpin C18-AX 100/pk
5010-21737	MonoSpin ME 50/pk
5010-21738	MonoSpin ME 100/pk
5010-01001	MonoSpin Custom Kit 20
5010-21740	MonoSpin Trial Kit1 10x4pk (C18, TiO, SCX, SAX)
5010-21741	MonoSpin Trial Kit2 10x4pk (C18, Amide, CBA, NH2)
5010-21742	MonoSpin Trial Kit3 10x4pk (SCX, SAX, CBA, NH2)